

RAB27A Antibody



PACO11643

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

RAB27A

Uniprot

P51159

Synonyms:

RAB27A, member RAS oncogene family; RAB27A; GS2; HsT18676; MGC117246; RAB27; RAM ;

Immunogen:

Human RAB27A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A