

RAB2B Antibody



PACO11646

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RAB2B

Uniprot

Q8WUD1

Synonyms:

RAB2B, member RAS oncogene family; RAB2B; FLJ14824 ;

Immunogen:

Human RAB2B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A