

RAB33A Antibody



PACO11649

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

RAB33A

Uniprot

Q14088

Synonyms:

RAB33A, member RAS oncogene family; RAB33A; MGC1488; RabS10 ;

Immunogen:

Human RAB33A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A