

RAB34 Antibody



PACO11650

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

RAB34

Uniprot

P0DI83

Synonyms:

RAB34, member RAS oncogene family; RAB34; RAB39; RAH ;

Immunogen:

Human RAB34.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A