

RAB3GAP1 Antibody



PACO11659

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

RAB3GAP1

Uniprot

Q15042

Synonyms:

RAB3 GTPase activating protein subunit 1 (catalytic); RAB3GAP1; DKFZp434A012; KIAA0066; P130; RAB3GAP; RAB3GAP130; WARBM1 ;

Immunogen:

Human RAB3GAP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A