

RABL2B Antibody



PACO11683

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RABL2B

Uniprot

Q9UNT1

Synonyms:

RAB, member of RAS oncogene family-like 2B; RABL2B; FLJ93981; FLJ98216 ;

Immunogen:

Human RABL2B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A