

RAD51AP1 Antibody



PACO11695

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

RAD51AP1

Uniprot

Q96B01

Synonyms:

RAD51 associated protein 1; RAD51AP1; PIR51 ;

Immunogen:

Human RAD51AP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A