

RAP1GDS1 Antibody



PACO11718

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RAP1GDS1

Uniprot

P52306

Synonyms:

RAP1, GTP-GDP dissociation stimulator 1; RAP1GDS1; GDS1; MGC118859; MGC118861; SmgGDS ;

Immunogen:

Human RAP1GDS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A