

RAPGEF5 Antibody



PACO11722

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RAPGEF5

Uniprot

Q92565

Synonyms:

Rap guanine nucleotide exchange factor (GEF) 5; RAPGEF5; GFR; KIAA0277; MR-GEF; REPAC ;

Immunogen:

Human RAPGEF5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A