

RAPGEFL1 Antibody



PACO11723

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

RAPGEFL1

Uniprot

Q9UHV5

Synonyms:

Rap guanine nucleotide exchange factor (GEF)-like 1; RAPGEFL1; Link-GEFII; MGC134798; MGC134799 ;

Immunogen:

Human RAPGEFL1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A