

RASSF8 Antibody



PACO11738

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

RASSF8

Uniprot

Q8NHQ8

Synonyms:

Ras association (RalGDS/AF-6) domain family (N-terminal) member 8; RASSF8; C12orf2; HoJ-1 ;

Immunogen:

Human RASSF8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A