

RBM12B Antibody



PACO11750

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RBM12B

Uniprot

Q8IXT5

Synonyms:

RNA binding motif protein 12B; RBM12B; MGC:33837 ;

Immunogen:

Human RBM12B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A