

# RBM19 Antibody



PACO11755

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RBM19

**Uniprot**

Q9Y4C8

**Synonyms:**

RNA binding motif protein 19; RBM19; DKFZp586F1023; KIAA0682 ;

**Immunogen:**

Human RBM19.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A