

RBM26 Antibody



PACO11756

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

RBM26

Uniprot

Q5T8P6

Synonyms:

RNA binding motif protein 26; RBM26; ARRS2; C13orf10; FLJ20957; MGC133295; MGC133296; PRO1777; RP11-255E21.1; SE70-2; ZC3H17 ;

Immunogen:

Human RBM26.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A