

REEP6 Antibody



PACO11797

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

REEP6

Uniprot

Q96HR9

Synonyms:

receptor accessory protein 6; REEP6; C19orf32; DP1L1; FLJ25383; TB2L1 ;

Immunogen:

Human REEP6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A