

# REG1B Antibody



PACO11798

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

REG1B

**Uniprot**

P48304

**Synonyms:**

regenerating islet-derived 1 beta; REG1B; PSPS2; REGH; REGI-BETA; REGL ;

**Immunogen:**

Human REG1B.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A