

RND2 Antibody



PACO11862

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

RND2

Uniprot

P52198

Synonyms:

Rho family GTPase 2; RND2; ARHN; RHO7; RhoN ;

Immunogen:

Human RND2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A