

RNF170 Antibody



PACO11869

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

RNF170

Uniprot

Q96K19

Synonyms:

ring finger protein 170; RNF170; DKFZp564A022; FLJ38306 ;

Immunogen:

Human RNF170.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A