

# RNMT Antibody



PACO11879

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RNMT

**Uniprot**

O43148

**Synonyms:**

RNA (guanine-7-) methyltransferase; RNMT; DKFZp686H1252; KIAA0398; MET; RG7MT1; hCMT1c ;

**Immunogen:**

Human RNMT.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A