

ROPN1L Antibody



PACO11888

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ROPN1L

Uniprot

Q96C74

Synonyms:

ropporin 1-like; ROPN1L; ASP; FLJ23003; FLJ25776 ;

Immunogen:

Human ROPN1L

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A