

RPL19 Antibody



PACO11914

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

RPL19

Uniprot

P84098

Synonyms:

ribosomal protein L19; RPL19; DKFZp779D216; FLJ27452; MGC71997 ;

Immunogen:

Human RPL19.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A