

RPL23 Antibody



PACO11916

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

RPL23

Uniprot

P62829

Synonyms:

ribosomal protein L23; RPL23; MGC111167; MGC117346; MGC72008; rpL17 ;

Immunogen:

Human RPL23.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A