

RPP38 Antibody



PACO11939

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

RPP38

Uniprot

P78345

Synonyms:

ribonuclease P/MRP 38kDa subunit; RPP38; RP11-455B2.5 ;

Immunogen:

Human RPP38.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A