

RPS10 Antibody



PACO11941

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RPS10

Uniprot

P46783

Synonyms:

ribosomal protein S10; RPS10; MGC88819 ;

Immunogen:

Human RPS10.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A