

RPS20 Antibody



PACO11948

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

RPS20

Uniprot

P60866

Synonyms:

ribosomal protein S20; RPS20; FLJ27451; MGC102930 ;

Immunogen:

Human RPS20.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A