

RPS27 Antibody



PACO11951

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

RPS27

Uniprot

P42677

Synonyms:

ribosomal protein S27; RPS27; MPS-1; MPS1 ;

Immunogen:

Human RPS27.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A