

RPS6KA5 Antibody



PACO11961

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

RPS6KA5

Uniprot

O75582

Synonyms:

ribosomal protein S6 kinase, 90kDa, polypeptide 5; RPS6KA5; MGC1911; MSK1; MSPK1; RLPK ;

Immunogen:

Human RPS6KA5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A