

RPS6KB2 Antibody



PACO11962

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

RPS6KB2

Uniprot

Q9UBS0

Synonyms:

ribosomal protein S6 kinase, 70kDa, polypeptide 2; RPS6KB2; KLS; P70-beta; P70-beta-1; P70-beta-2; S6K-beta2; S6K2; SRK; STK14B; p70(S6K)-beta; p70S6Kb ;

Immunogen:

Human RPS6KB2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A