

RSPO2 Antibody



PACO11975

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RSPO2

Uniprot

Q6UXX9

Synonyms:

R-spondin 2 homolog (Xenopus laevis); RSPO2; CRISTIN2; MGC35555; MGC43342 ;

Immunogen:

Human RSPO2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A