

RTP1 Antibody



PACO11984

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RTP1

Uniprot

P59025

Synonyms:

receptor (chemosensory) transporter protein 1; RTP1; MGC35450 ;

Immunogen:

Human RTP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A