

# SALL2 Antibody



PACO12018

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

SALL2

**Uniprot**

Q9Y467

**Synonyms:**

sal-like 2 (Drosophila); SALL2; FLJ10414; FLJ55746; HSAL2; KIAA0360; ZNF795; p150(Sal2) ;

**Immunogen:**

Human SALL2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A