

# SAMM50 Antibody



PACO12022

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SAMM50

**Uniprot**

Q9Y512

**Synonyms:**

sorting and assembly machinery component 50 homolog (S. cerevisiae); SAMM50; CGI-51; FLJ35825; FLJ42905; FLJ99036; OMP85; SAM50; TOB55; TRG-3; YNL026W ;

**Immunogen:**

Human SAMM50.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A