

SAP18 Antibody



PACO12024

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

SAP18

Uniprot

O00422

Synonyms:

Sin3A-associated protein, 18kDa; SAP18;2HOR0202; MGC27131; SAP18P ;

Immunogen:

Human SAP18.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A