

SCARF1 Antibody



PACO12042

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SCARF1

Uniprot

Q14162

Synonyms:

scavenger receptor class F, member 1; SCARF1; KIAA0149; MGC47738; SREC ;

Immunogen:

Human SCARF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A