

SCFD2 Antibody



PACO12045

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SCFD2

Uniprot

Q8WU76

Synonyms:

sec1 family domain containing 2; SCFD2; FLJ21060; FLJ39514; STXBP1L1 ;

Immunogen:

Human SCFD2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A