

SCUBE3 Antibody



PACO12062

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SCUBE3

Uniprot

Q8IX30

Synonyms:

signal peptide, CUB domain, EGF-like 3; SCUBE3; CEGF3; DKFZp686B09105; DKFZp686B1223; DKFZp686D20108; FLJ34743 ;

Immunogen:

Human SCUBE3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A