

SDCBP2 Antibody



PACO12067

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SDCBP2

Uniprot

Q9H190

Synonyms:

syndecan binding protein (syntenin) 2; SDCBP2; FLJ12256; SITAC18; ST-2 ;

Immunogen:

Human SDCBP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A