

NEMF Antibody



PACO12068

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

NEMF

Uniprot

O60524

Synonyms:

serologically defined colon cancer antigen 1; SDCCAG1; FLJ10051; NY-CO-1 ;

Immunogen:

Human SDCCAG1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A