

CWC27 Antibody



PACO12069

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CWC27

Uniprot

Q6UX04

Synonyms:

serologically defined colon cancer antigen 10; SDCCAG10; NY-CO-10 ;

Immunogen:

Human SDCCAG10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A