

SDCCAG3 Antibody



PACO12070

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SDCCAG3

Uniprot

Q96C92

Synonyms:

serologically defined colon cancer antigen 3; SDCCAG3; NY-CO-3 ;

Immunogen:

Human SDCCAG3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A