

SEC23IP Antibody



PACO12082

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SEC23IP

Uniprot

Q9Y6Y8

Synonyms:

SEC23 interacting protein; SEC23IP; MSTP053; P125 ;

Immunogen:

Human SEC23IP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A