

SEC24A Antibody



PACO12083

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

SEC24A

Uniprot

O95486

Synonyms:

SEC24 family, member A (S. cerevisiae);

Immunogen:

Human SEC24A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A