

SEC24C Antibody



PACO12084

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SEC24C

Uniprot

P53992

Synonyms:

SEC24 family, member C (S. cerevisiae); SEC24C; KIAA0079 ;

Immunogen:

Human SEC24C.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A