

SEMA4B Antibody



PACO12094

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

SEMA4B

Uniprot

Q9NPR2

Synonyms:

sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 4B; SEMA4B; KIAA1745; MGC131831; SEMAC; SemC ;

Immunogen:

Human SEMA4B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A