

SERINC3 Antibody



PACO12103

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SERINC3

Uniprot

Q13530

Synonyms:

serine incorporator 3; SERINC3; AIGP1; DIFF33; SBB199; TDE; TDE1; TMS-1 ;

Immunogen:

Human SERINC3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A