

SERPINA10 Antibody



PACO12106

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SERPINA10

Uniprot

Q9UK55

Synonyms:

serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 10; SERPINA10; PZI; ZPI ; Serpin A10

Immunogen:

Human SERPINA10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A