

SERPINA5 Antibody



PACO12108

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SERPINA5

Uniprot

P05154

Synonyms:

serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 5; SERPINA5; PAI3; PCI; PLANH3; PROCI ;

Immunogen:

Human SERPINA5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A