

SERPINB12 Antibody



PACO12111

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SERPINB12

Uniprot

Q96P63

Synonyms:

serpin peptidase inhibitor, clade B (ovalbumin), member 12; SERPINB12; MGC119247; MGC119248; YUKOPIN ;

Immunogen:

Human SERPINB12.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A