

SUGP1 Antibody



PACO12139

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SUGP1

Uniprot

Q8IWZ8

Synonyms:

splicing factor 4; SF4; DKFZp434E2216; F23858; RBP ;

Immunogen:

Human SF4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A