

AKAP17A Antibody



PACO12146

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

AKAP17A

Uniprot

Q02040

Synonyms:

splicing factor, arginine/serine-rich 17A; SFRS17A;721P; CCDC133; CXYorf3; DXYS155E; MGC125365; MGC125366; MGC39904; XE7; XE7Y ;

Immunogen:

Human SFRS17A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A